



## 8. Analysis of a Mixture of Carbonate and Bicarbonate




Green Profile  
See Section 0


This procedure involves two titrations. First, total alkalinity ( $[HCO_3^-] + 2[CO_3^{2-}]$ ) is measured by titrating the mixture with standard HCl to a bromocresol green end point:



A separate aliquot of unknown is treated with excess standard NaOH to convert  $HCO_3^-$  to  $CO_3^{2-}$ :



Then all the carbonate is precipitated with  $BaCl_2$ :



The excess NaOH is immediately titrated with standard HCl to determine how much  $HCO_3^-$  was present. From the total alkalinity and bicarbonate concentration, you can calculate the original carbonate concentration.

### Reagents

*Standard 0.1 M NaOH and standard 0.1 M HCl:* From Experiment 6.

*CO<sub>2</sub>-free water:* Boil 500 mL of distilled water to expel CO<sub>2</sub> and pour the water into a 500-mL plastic bottle. Screw the cap on tightly and allow the water to cool to room temperature.

Keep tightly capped when not in use.

*10 wt% aqueous BaCl<sub>2</sub>:* 35 mL/student.

*Bromocresol green and phenolphthalein indicators:* See Experiment 7 for recipes.

*Unknowns:* Solid unknowns (2.5 g/student) can be prepared from reagent-grade sodium carbonate or potassium carbonate and bicarbonate. Unknowns should be stored in a desiccator and should not be heated. Heating at 50°–100°C converts NaHCO<sub>3</sub> to Na<sub>2</sub>CO<sub>3</sub>.

**Procedure**

1. Accurately weigh 2.0–2.5 g of unknown into a 250-mL volumetric flask by weighing the sample in a capped weighing bottle, delivering some to a funnel in the volumetric flask, and reweighing the bottle. Continue this process until the desired mass of reagent has been transferred to the funnel. Rinse the funnel repeatedly with small portions of CO<sub>2</sub>-free water to dissolve the sample. Remove the funnel, dilute to the mark, and mix well.
2. *Total alkalinity:* Pipet a 25.00-mL aliquot of unknown solution into a 250-mL flask and titrate with standard 0.1 M HCl, using bromocresol green indicator as in Experiment 6 for standardizing HCl. Repeat this procedure with two more 25.00-mL aliquots.
3. *Bicarbonate content:* Pipet 25.00 mL of unknown and 50.00 mL of standard 0.1 M NaOH into a 250-mL flask. Swirl and add 10 mL of 10 wt% BaCl<sub>2</sub>, using a graduated cylinder. Swirl again to precipitate BaCO<sub>3</sub>, add 2 drops of phenolphthalein indicator, and immediately titrate with standard 0.1 M HCl. Repeat this procedure with two more 25.00-mL samples of unknown.
4. From the results of step 2, calculate the total alkalinity and its standard deviation. From the results of step 3, calculate the bicarbonate concentration and its standard deviation. Using the standard deviations as estimates of uncertainty, calculate the concentration (and uncertainty) of carbonate in the sample. Express the composition of the solid unknown in a form such as 63.4 (±0.5) wt% K<sub>2</sub>CO<sub>3</sub> and 36.6 (±0.2) wt% NaHCO<sub>3</sub>.